Product Data Sheet

ESCAPIES Evacuation mat Page 1 of 1

GENERAL INFORMATION

REF	PRODUCT DESCRIPTION	C = MP Class R = Rule	UMDNS	РНОТО
5000	The evaculation mat that is placed between the mattress and the bed slats is allowing nursing home residents to be evacuated by a single staff member. Bed-ridden residents lying on an evacuation mat can be slid across a variety of surfaces and through narrow staircases.	CI - R1	15831	

Main Fabric:	100% Polyester
Dyes:	Pigment & Disperse Dyes
Drawstrings:	PP (Polypropylene)
Elastic:	100% Polyester & Rubber
Buckle:	PP (Polypropylene)
Coating:	Polyurethane Coating
Flame Retardant:	ISO 6940
Load Capacity:	150 kg
Dimensions:	195 x 93 cm

Headquarters: Dutchblue.com LTD

26 Parizhka Komuna str. E info@dutchblue.com W www.dutchblue.com 9000 Varna, Bulgaria

